Seraclone™ Anti-A (AB01), -B (AB02), -AB (AB03)

For typing of AB0 blood characteristics

REF	801 320	10 ml	Seraclone™ Anti-A
REF	801 325	10x10 ml	Seraclone™ Anti-A
REF	801 345	10 ml	Seraclone™ Anti-B
REF	801 350	10x10 ml	Seraclone™ Anti-B
REF	801 370	10 ml	Seraclone™ Anti-AB
REF	801 375	10x10 ml	Seraclone™ Anti-AB
PRES	Preservative: < 0.1 % NaN ₃		
ACT	Activity / Min.titer		

IVD Diagnostic reagent for in vitro use only

To be used by trained laboratory personnel only. The **ACT** is printed on the label. Symbols are used according to regulation EN 980. Symbols not contained in regulation EN 980 are described.

Test purpose

AB0 blood type characteristics are determined with the monoclonal reagents Anti-A, Anti-B and Anti-AB. They are verified by serum characteristics determination. Blood typing is considered valid if both red blood cell and serum characteristics were examined.

Test principle

The test principle is a hemagglutination test for the detection of red blood cell antigens. The antibodies in Seraclone™ Anti-A, -B, -AB bind to the corresponding antigen on red blood cells and cause an antigen-antibody-reaction visible as red blood cell agglutination. The four AB0 blood types A, B, AB and 0 are defined by the presence or absence of A and B characteristics on red blood cells. The absence of both A and B characteristics defines blood type 0. The antigens (characteristics) A and B react with the corresponding antibody in Seraclone™ Anti-A, -B, -AB. Seraclone™ Anti-A, -B, -AB can be used for slide, tile or tube testing.

Reagent

Seraclone™ Anti-A, -B, -AB contain as reactive components monoclonal antibodies of the immunoglobulin class IgM. They are derived from hybridoma cell lines which are created by fusing mouse antibody producing B-lymphocytes with mouse myeloma cells and demonstrate consistent specificity and reproducibility characteristic for monoclonal antibodies. Both antibodies derived from a single clone (sister cells of one hybridoma cell) and a mixture of different antibodies derived from several clones are called monoclonal. The bovine albumin used for the production of this reagent is purchased from BSE-free sources.

Seraclone™ Anti-A	Clone A003
Seraclone™ Anti-B	Clone B005
Seraclone™ Anti-AB	Clones BS 63/BS 85

Material required but not supplied

- Slides, Bioplate
- Pipettes (drop volume 40-50µl)
- Isotonic saline or Phosphate Buffered Solution (PBS). (PBS only when tested with Seraclone™ Anti-AB)
- Glass tubes
- Laboratory centrifuge
- Seraclone Control AB0+Rh (REF 805171)
- Control Set QC (REF 816034)

Specimen collection

Blood samples should be taken following general blood sampling guidelines. Fresh, non-hemolytic samples (native, EDTA or citrate blood) should be used. The specimen should be tested as soon as possible after collection. If testing is delayed, specimens should be stored at 2...8°C. If needed prepare a red cell suspension from a centrifugated blood sample (e. g. 2 min. at 1000 x g) in isotonic saline as suspension medium. ◀ Strongly lipemic, icteric or microbiologically contaminated samples may lead to false results.

Test procedure

For the detection of AB0 blood types the red cells to be tested are treated with Seraclone™ Anti-A, Anti-B, Anti-AB. The serum to be tested is treated with reagent red cells, such as Biotestcell™-A₁, -A₂, -B, -0 (reverse typing).

Rapid test

- Place 1 drop reagent to slide (do not pre-heat!).
- Add 1 small drop whole blood and mix well.
- Agglutination will occur within 30 to 60 seconds. In order not to overlook weak A reactions, the agglutination should be read only after 2 minutes of careful rotation.

Incubation test

- Suspend red cells to be tested 5-10% in isotonic saline solution or use whole blood.
- On a slide mix 1 drop reagent and 1 drop red cell suspension or whole blood.
- Incubate 15-30 minutes at room temperature (cover adequately to prevent drying out)
- Gently rock slide and observe for agglutination.

Tube test

- Suspend red cells to be tested in 3-5% isotonic saline or PBS*.
- In a properly marked tube mix 1 drop reagent and 1 drop red blood cell suspension.
- Centrifugate 2 minutes at 150 - 200 x g or 20 seconds at 800 - 1000 x g or incubate 20 minutes at room temperature.
- Gently dislodge cell button and observe for agglutination.

*PBS is only approved for Seraclone™ Anti-AB

Test of serum characteristics (reverse typing)

The serum to be tested should be examined with known A₁-, A₂-, B- and 0-cells or Biotestcell™-A₁, -A₂, -B, -0 with the above described procedures in slide, tile or tube testing. Since serum characteristics may react in different strength, an incubation for 15-30 minutes at room temperature should be performed. Generally, newborns and young babies do not show test reaction due to missing isoagglutinines. Isoagglutinines may also be absent in elderly patients.

Interpretation of results

Reaction patterns red cell characteristics and isoagglutinines

Test sera + Patient red cells			Test red cells + serum/plasma				Blood- group
Anti-A	Anti-B	Anti-AB	A ₁	A ₂	B	0	
+	-	+	-	-	+	-	A
-	+	+	+	+	-	-	B
-	-	-	+	+	+	-	0
+	+	+	-	-	-	-	AB

+ = agglutination - = no agglutination

Reactions of Seraclone™ Anti-A, Anti-B and Anti-AB with AB0 variations

Seraclone™ test reagents do not react with cryptoantigens (T-, Tn-, Tk-activated cells).

Seraclone™ Anti-B reacts correctly negative with acquired B characteristics.

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Evaluation of the reaction strength is carried out according to the Technical Manual (1):

Reaction strength	Agglutination
4+	One single agglutinate
3+	Several large agglutinates
2+	Medium-size agglutinates, clear background
1+	Small agglutinates, turbid background
+/-	Barely visible agglutination, turbid background
-	No agglutination

Quality control

The reactivity of all blood typing reagents should be confirmed on each day of use. To confirm the reactivity or specificity of Bio-Rad monoclonal AB0 reagents (Anti-A, Anti-B and/or Anti-AB) each should be tested with antigen-positive and antigen-negative red blood cells, respectively. Each reagent is satisfactory for use if it reacts only with antigen-positive red blood cells. Confirmation of results in forward grouping must be obtained by performing the reverse grouping test.

Performance characteristics and limitations of the method

In very rare cases, it has been observed that when using Seraclone™ Anti-A minor traces of A-blood group substance were detected on B-erythrocytes. If this B (A) phenomenon does occur, it should be solved like other discrepancies between blood group determination and reverse typing.

Seraclone™ Anti-B reacts correctly negative with acquired B characteristics.

In case of ambiguous results it is recommended to wash red blood cells at least 2 times.

In case of questionable results of unknown origin our Bio-Rad Service will assist you (phone: +49-6103-3130-611).

Shelf life

Storage temperature of the Bio-Rad monoclonal AB0 reagents is 2...8°C. After opening the reagent, the product can be stored until its expiry date under proper storage conditions (2...8°C). Recurrent alternating storage of the

reagent Seraclone™ Anti-AB from 2...8°C to ≥ room temperature may result in protein precipitation, which does not affect the efficacy of the reagent. The expiry date is printed on the label.

Warning and precautions

- Do not use damaged vials.
- Turbidity or other visible changes may indicate a bacterial contamination. Do not use if turbid.
- Manual techniques are to be performed according to the manufacturer's instructions. Each deviation from these instructions is the sole responsibility of the user.
- Seraclone™ Anti-A, -B, -AB is suited for use in automated systems. The use in automated systems is to be validated by the user. The use of Seraclone™ Anti-A, -B, -AB in automated systems may require dilution of the reagent. The user's responsibility extends to all changes in the reagent, e.g. freezing on microtiterplates.
- Used tests must be discarded as hazardous material. Waste management information can be found in the safety data sheet.
- Because of biotechnological manufacturing of non-human monoclonal antibodies, an infection with hepatitis, HIV 1/2 or other infectious diseases can be considered impossible. Nevertheless, all test reagents of biological origin must be regarded as potential transmitters of hepatitis, HIV or other infectious agents. Appropriate safety precautions are recommended.
- The packaging of this product contains natural rubber latex which may cause allergic reactions.
- Seraclone™ Anti-B contains tartrazine which may cause allergic reaction.
- The products contain sodium azide (NaN₃), which may react with lead or copper plumbing to form explosive azides. If discarded in the sink, flush with large amounts of water to prevent the build-up of explosive metal azides.

References

(1) Technical Manual, 17th edition, Section 1, American Association of Blood Banks

Key: Underline = Addition or significant change ◀ = Deletion of text

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Reazioni di Seraclone™ Anti-A, Anti-B e Anti-AB con varianti AB0